

**General lysate preparation**

- ↓ 50 mM NaOH (600 µl for ~ 2 mm tail) (Adjust the volume of NaOH if you have a smaller tissue)
- ↓ Dissolved the tail at 100°C for 20-30 min
  - (Use cap to lock the eppendorf. Invert or shake the eppendorf sometimes)
- ↓ Spin at top speed for 1 min.
- ↓ Neutralized with 50 µl of Tris-HCl (1M, pH 8.0). Mix well (or vortex) and Spin again.
- ↓ put on ice and ready for used for genotyping.

**From Web 1 (<http://research.uci.edu/tmf/DNAprep.htm>)**

Final concentration of tail digestion buffer (TDB):

<b>Chemical</b>	<b>Final Conc</b>	<b>200 ml</b>	
KCl	50 mM	0.7455 g	
Tris-HCl (pH 9.0)	10 mM	0.3152 g	Using 1N NaOH adjust pH
Triton X-100	0.1 %	0.2 ml	

\* Before use, add Proteinase K (stock 20 mg/ml) 5 µl /ml

\*\* Proteinase K from Novagen (cat#70663-4), prepare with 50 mM Tris-HCl, pH8.0 with 1 mM CaCl<sub>2</sub>.

**Procedure:**

- ↓ Add 100 µl of TDB per tail piece (2-5 mm) in a microcentrifuge tube. Be careful not to cut too much tail.
- ↓ Incubate the tube at either 60°C for 3 hours with gentle mixing or 55°C overnight.
- ↓ Incubate the tube at 95°C (or boil) for 10 minutes to denature the Proteinase K.
- ↓ Spin in microcentrifuge at top speed for 15 minutes. Use an aliquot of supernatant straight from the tube (e.g., 0.5 µl in a 10 µl reaction) for PCR. If the results are questionable, try a 1:20 dilution of the DNA

**Genotyping for L4-CRE line (JAX no.009613);**

**dlx6a-CRE line (JAX no. 008199); ROR-CRE line**

**Primer:** oIMR1084: 5'-GCGGTCTGGCAGTAAAAACTATC-3'  
oIMR1085: 5 '-GTGAAACAGCATTGCTGTCACCT-3'  
oIMR7338 (control): 5'- CTAGGCCACAGAATTGAAAGATCT-3'  
oIMR7339 (control): 5'- GTAGGTGGAAATTCTAGCATCATCC-3'

\* Primers are from Jackson Lab Web

**Reaction:**

	Volume (μl)
2 x Econo master mix	5
Primer-oIMR1084 (10 μM)	1
Primer-oIMR1085(10 μM)	1
Primer oIMR7338 (control)	1
Primer oIMR7339 (control)	1
Template (Genomic DNA)	1
Final	10

**Use filter pipette tips when aliquoting reagents for PCR as well as when preparing master mix!**

**PCR program:**      95°C, 5 min  
                        95°C, 30 sec / 52°C, 40 sec / 72°C, 40 sec (35 cycles)  
                        72°C, 5 min

**Gel:** 3% agarose gel (**necessary**), 135, 40 min

**Results:** This assay will NOT distinguish hemizygous from homozygous transgenic animals.

Expected results: Transgene = ~100 bp; internal positive control= 324 bp  
and if see transgene express, the internal positive control will be weak

**Protocol-Genotyping**  
**Jessica Felker/Dr. Hui-Chen Lu's Lab**

**Genotyping of CB1 floxed/floxed**

**Primer:** G50/CB50 : 5'-CTGTCTCTGGTCCTCTTAAA-3'  
G51/CB51: 5'-GGTGTACACCTCTGAAACAGA-3'

**Reaction:**

	Volume ( $\mu$ l)
H <sub>2</sub> O	2
2 x Econo master mix	5
Primer-G50 (10 $\mu$ M)	1
Primer-G51 (10 $\mu$ M)	1
Template (Genomic DNA)	1
Final	10

**Use filter pipette tips when aliquoting reagents for PCR as well as when preparing master mix!**

**PCR program:**      95°C, 5 min  
                        95°C, 30 sec / 58°C, 40 sec / 72°C, 40 sec (35 cycles)  
                        72°C, 5 min

**Gel:** 2% agarose gel, 135 45 min

**Results:** CB<sub>1</sub> locus:    +/+ : 400 bp  
                              f/+ : 500 bp and 400 bp  
                              f/f : 500 bp

**Protocol-Genotyping**  
**Jessica Felker/Dr. Hui-Chen Lu's Lab**

**Genotyping of CB1 knockout**

**Primer:** CB50: 5'-GCTGTCTCTGGTCCTCTTAAA-3'  
CB51: 5'-GGTGTACACCTCTGAAAACAGA-3'  
CB54: 5'-CCTACCCGGTAGAATTAGCTT-3'

CB-WT primer mix: CB50 + CB51

CB-neo primer mix: CB50 + CB54

**Reaction:**

	Volume ( $\mu$ l)
H <sub>2</sub> O	2
2 x Econo master mix	5
Primer-CB50 (10 $\mu$ M)	1
Primer-CB51(for WR) or CB54 (for neo) (10 $\mu$ M)	1
Template (Genomic DNA)	1
Final	10

**Use filter pipette tips when aliquoting reagents for PCR as well as when preparing master mix!**

**PCR program:** 95°C, 5 min  
95°C, 30 sec / 58°C, 40 sec / 72°C, 40 sec (35 cycles)  
72°C, 5 min

**Gel:** 2% agarose gel, 135V, 45 min

**Results:**

+/+ : 400 bp

**Genotyping of VGAT-CRE line; ROR $\alpha$ -CRE line; NEX-Cre line;**  
**NEX-Cre-ER line**

**Primer:** PCR Cre484 : 5'- GCATTTCTGGGGATTGCTTA -3'  
PCR Cre834: 5'-GTCATCCTAGCGCCGTAAA -3'

**Reaction:**

	Volume ( $\mu$ l)
H <sub>2</sub> O	2
2 x Econo master mix	5
Primer- Cre484(10 $\mu$ M)	1
Primer- Cre834(10 $\mu$ M)	1
Template (Genomic DNA)	1
Final	10

**Use filter pipette tips when aliquoting reagents for PCR as well as when preparing master mix!**

**PCR program:** 95°C, 5 min  
95°C, 30 sec / 62°C, 30 sec / 72°C, 20 sec (35 cycles)  
72°C, 5 min

**Gel:** 2% agarose gel, 135V, 30 min

**Results:** the primer were designed in the coding sequence of CRE,  
the amplicon size is 351 bp

**Protocol-Genotyping**  
**Jessica Felker/Dr. Hui-Chen Lu's Lab**

**Genotyping of NEX-CRE line; NEX-CRE-ERT2 line**

NEX-Cre = knock-in of Cre ORF into the NEX gene locus

\*\*Note that in the homozygous state, this effectively knocks-out NEX function\*\*

**Primer:**

NexForward: 5'- GAG TCC TGG AAT CAG TCT TTT TC -3'


Localized 148bp in front of the start ATG

NexReverse: 5'-AGA ATG TGG AGT AGG GTG AC-3'

Localized 622bp behind the start ATG

CreReverse: 5'-CCG CAT AAC CAG TGA AAC AG-3'

Localized 370bp behind the start ATG, within the Cre coding region



**Reaction:**

	Volume (μl)
H <sub>2</sub> O	2
2 x Econo master mix	5
NEX-forward (10 μM)	1
NEX-reverse (for wild type)/ CRE-reverse (for CRE) (10 μM)	1
Template (Genomic DNA)	1
Final	10

**Use filter pipette tips when aliquoting reagents for PCR as well as when preparing master mix!**

**PCR program:** 95°C, 5 min

95°C, 30 sec / 54°C, 30 sec / 72°C, 60 sec (40 cycles)

72°C, 5 min

**Gel:** 1.5 % agarose gel, 135V, 40 min

**Results:**

Wildtype allele = 770bp product of NexForward and NexReverse primers

NEX-Cre allele = 520bp product of NexForward and CreReverse primers

**Protocol-Genotyping**  
Jessica Felker/Dr. Hui-Chen Lu's Lab

**Genotyping of mGluR5 flox (MGF) line**

**Primer:** Forward Primer 153: 5'-AGATGTCCCACCTACCTGATGT-3'

Reverse Primer 154: 5'- AGTTCCGTGTCTTATTCTTAGC-3'

**Reaction:**

	Volume ( $\mu$ l)
H <sub>2</sub> O	2
2 x Econo master mix	5
Primer-MGF-F (10 $\mu$ M)	1
Primer-MGF-R(10 $\mu$ M)	1
Template (Genomic DNA)	1
Final	10

**Use filter pipette tips when aliquoting reagents for PCR as well as when preparing master mix!**

**PCR program:** 95°C, 5 min

95°C, 30 sec / 52°C, 30 sec / 72°C, 30 sec (35 cycles)

72°C, 5 min


**Gel:** 3% agarose gel, 135V, 40 min

**Results:**

Image (2).JPG

**Result:**

wt: ~200 bp (primer 153~154)  
flox: ~250 bp (primer 153~154)  
deletion: ~360 bp (primer 153~316)



**Protocol-Genotyping**  
Jessica Felker/Dr. Hui-Chen Lu's Lab

**Genotyping of mGluR5 knockout line (JAX no. 003558)**

**Primer:** oIMR1034 (WT-F): 5'- CAC ATg CCA ggT gAC ATC AT -3'  
oIMR1035 (WT-R): 5'- CCA TgC Tgg TTg CAg AgT AA -3'  
oIMR0013 (neo-F): 5'- CTT ggg Tgg AgA ggC TAT TC -3'  
oIMR0014 (neo-R): 5'- Agg TgA gAT gAC Agg AgA TC -3'

**Reaction:**

	Volume (μl)
H <sub>2</sub> O	2.5
2 × Econo master mix	5
Primer-F (10 μM)	1
Primer-R(10 μM)	1
Template (Genomic DNA)	0.5
Final	10

**Use filter pipette tips when aliquoting reagents for PCR as well as when preparing master mix!**

**PCR program:** 95°C, 5 min  
95°C, 30 sec / 59°C, 25 sec / 72°C, 30 sec (30 cycles)  
72°C, 5 min

**Gel:** 1.5% agarose gel, 135V, 40 min

**Results:** wt band 442 bp and neo band 280 bp

**Genotyping of tom/tom mice**

stock number of Jackson Lab: 007909

Allele: Gt (ROSA)26Sor tm9<sup>(CAG-tdTomato)Hze</sup>

**Primer:** oIMR9020 (WT-F): 5'-AAG GGA GCT GCA GTG GAG TA-3'

oIMR9021 (WT-R): 5'-CCG AAA ATC TGT GGG AAG TC-3'

oIMR9103 (Mutant-R): 5'-GGC ATT AAA GCA GCG TAT CC-3'

oIMR9105 (Mutant-F): 5'-CTG TTC CTG TAC GGC ATG G-3'

**Reaction:**

	Volume ( $\mu$ l)
H <sub>2</sub> O	0.5
2 × Econo master mix	5
Primer-oIMR9020 (10 $\mu$ M)	1
Primer-oIMR9021(10 $\mu$ M)	1
Primer oIMR9103 (control)	1
P Primer oIMR9105 (control)	1
Template (Genomic DNA)	0.5
Final	10

**Use filter pipette tips when aliquoting reagents for PCR as well as when preparing master mix!**

**PCR program:** 95°C, 5 min

95°C, 30 sec / 61°C, 30 sec / 72°C, 25 sec (35 cycles)

72°C, 5 min

**Gel:** 2% agarose gel, 135V, 40min, should run longer to separate

**Results:** Mutant = 196 bp; Heterozygote = 297 bp and 196 bp; Wild Type = 297 bp

**Protocol-Genotyping**  
**Jessica Felker/Dr. Hui-Chen Lu's Lab**

**Genotyping of tau-stop-mGFP mice**

**Primer:**

GFP-F: CGGCGAGGGCGAGGGCGATG

GFP-R: CAGGGGGCCGTCGCCGATGG

It expresses a membrane bound form of eGFP (N-terminal 40Aa of MARCKS) and is followed by an IRES-NLS-LacZ. The membrane bound version of eGFP is very good to see all details of cell surface.

Tau<sup>mGFP</sup> mouse is a Cre-reporter line containing a floxed ‘stop transcription’ sequence in front of membrane-anchored green fluorescent protein (mGFP) and an IRES-NLS-lacZ gene inserted into exon 2 of the Tau locus ([Hippenmeyer et al., 2005](#)).

Cre-mediated recombination can be detected by the presence of nuclear β-galactosidase, because of the nuclear localization sequence (NLZ) engineered in the LacZ gene, and by the expression of mGFP in the axons of recombined neurons.

*Hippenmeyer, S., Vrieseling, E., Sigrist, M., Portmann, T., Laengle, C., Ladle, D.R. & Arber, S. (2005) A developmental switch in the response of DRG neurons to ETS transcription factor signaling. PLoS Biol., 3, e159.*

**Reaction:**

	Volume (μl)
H <sub>2</sub> O	2.5
2 × Econo master mix	5
Primer-GFP-F (10 μM)	1
Primer- GFP-R (10 μM)	1
Template (Genomic DNA)	0.5
Final	10

**Use filter pipette tips when aliquoting reagents for PCR as well as when preparing master mix!**

**PCR program:**      95°C, 5 min  
                          95°C, 30 sec / 65°C, 40 sec / 72°C, 40 sec (30 cycles)  
                          72°C, 5 min

**Gel:** 1.5% agarose gel, 135V, 40 min

**Results:** 500 bp transgene; no band is wild type

ROR $\alpha$ -Cre mice were generated by inserting an *IRES-cre* cDNA fragment into the 3' noncoding region of the ROR $\alpha$  gene in Dr Dennis O'Leary's laboratory (data not shown). Cre expression in ROR $\alpha$ -Cre mice is similar to endogenous ROR $\alpha$  expression ([Nakagawa & O'Leary, 2003](#)).

*Nakagawa, Y. & O'Leary, D.D. (2003) Dynamic patterned expression of orphan nuclear receptor genes ROR $\alpha$  and ROR $\beta$  in developing mouse forebrain. Dev. Neurosci., 25, 234–244.*

**Genotyping of NMNAT2 mice**

**Primer:**

The integration site is on chr1G3 and is linked to the right IR/DR of the transposon beginning at 154,851,334(+)

**R3 primer** (from the right IR/DR): 5' - CCACTGGGA ATG TGA TGA AAG AAA TAA AAG C -3'  
(ove2172Chr1G3)

**RF (right flank) primer.** This primer is antisense and is downstream from the integration site  
5' -CTG ACG TCT ATC TAG AAG TAC AC -3' (chr1G3:154,851,749-154,851,771(-))

ove2172Chr1G3 **upstream primer (sense)**

5'- GGA AAA GCC CAA ATCCCA GCA A -3' (chr1G3:154,850,990-154,851,011(+))

**Reaction**

	Volume ( $\mu$ l)
H <sub>2</sub> O	2
2 x Econo master mix	5
RF (10 $\mu$ M)	1
R3 or UP (10 $\mu$ M)	1
Template (Genomic DNA)	1
Final	10

**Use filter pipette tips when aliquoting reagents for PCR as well as when preparing master mix!**

**PCR program:**      95°C, 5 min  
                          95°C, 30 sec / 60°C, 40 sec / 72°C, 40 sec (35 cycles)  
                          72°C, 5 min

**Gel:** 2 % agarose gel, 135V, 40h

**Results:**

R3 + RF: 630 bp band if the mouse is transgenic; no band for a non-transgenic mouse.  
RF + upstream primer: 760 bp band if the mouse is wild type or heterozygous; no band if the mouse is homozygous.

**Male-specific Sry (Sex determining region protein gene) gene**

**Primer:**

Sry-F (Forward locus 8276-8295): 5'- TGGGACTGGTGACAATTGTC -3'

Sry-R (Reverse locus 8677-8658): 5'- GAGTACAGGTGTGCAGCTCT -3'

\*from Journal of Neuroscience Methods 95 (2000), 127-132. Quick sex determination of mouse fetuses

**Reaction: use Econo 2x master mix**

	Volume ( $\mu$ l)
H <sub>2</sub> O	2.5
2 × master mix	5
Sry-F (10 $\mu$ M)	1
Sry-R (10 $\mu$ M)	1
Template (Genomic DNA)	0.5
Final	10

**Use filter pipette tips when aliquoting reagents for PCR as well as when preparing master mix!**


**PCR program:**      95°C, 5 min  
                        95°C, 30 sec / 57°C, 40 sec / 72°C, 40 sec (35 cycles)  
                        72°C, 5 min

**Gel:** 2 % agarose gel, 135V, 40h

**Results:** 402 bp representing the Y-specific Sry gene

### Genotyping of FGF9 conditional KO mice

**Mouse is from Fen Wang (MGI ID: 3621452)**



#### **Primer:**

FGF9-P3 (Forward): AGA GAA ACT GCC CTG TCC AAC CAA

FGF9-P4 (Reverse) : CAG CCC GAA GAC ATT CGG CCA CAA

#### **Reaction**

	Volume ( $\mu$ l)
H <sub>2</sub> O	2.5
Econo 2x mater mix	5
Primer- F (10 $\mu$ M)	1
Primer- R(10 $\mu$ M)	1
Template (Genomic DNA)	0.5
Final	10

**Use filter pipette tips when aliquoting reagents for PCR as well as when preparing master mix!**


**PCR program:**      95°C, 5 min  
 95°C, 30 sec / 60°C, 30 sec / 72°C, 30 sec (34 cycles)  
 72°C, 5 min

**Gel:** 2 % agarose gel, 135V, 30 min

**Results:** WT allel 170 bp; flox allel 383bp

### Genotyping of FGFR2 conditional KO mice

Original paper Development, 2003, 130, 3063-3074



#### **Primer:**

FGFR2-F1 (Forward): ATAGGAGAACAGGGCGG

FGFR2-R1 (Reverse) : TGCAAGAGGCACCAGTCAG

#### **Reaction**

	Volume (μl)
H <sub>2</sub> O	2.5
Econo 2x mater mix	5
Primer- F (10 μM)	1
Primer- R(10 μM)	1
Template (Genomic DNA)	0.5
Final	10

**Use filter pipette tips when aliquoting reagents for PCR as well as when preparing master mix!**

**PCR program:** 95°C, 5 min

95°C, 30 sec / 60°C, 30 sec / 72°C, 30 sec (34 cycles)

72°C, 5 min

\* the same program as FGF9

**Gel:** 2 % agarose gel, 135V, 30 min

**Results:** WT allel 142 bp; flox allel 207 bp

**Protocol-Genotyping**  
**Jessica Felker/Dr. Hui-Chen Lu's Lab**

**Genotyping of Ntrk1 mice (JAX no. 022362)**

These mice carry a hypomorphic floxed F592A mutation in mouse *Ntrk1* (neurotrophic tyrosine kinase, receptor, type 1; also called TrkA). Signaling is effectively blocked by application of 1NMPP1. This strain is useful in studies of nerve growth and survival.

Original paper: Neurons, 2005, Volume 46, Issue 1, pages 13-21.  
A Chemical-Genetic Approach to Studying Neurotrophin Signaling

**Primer:**

Ntrk1-F (16834): 5'-CAC AGG GGC TGG AAA CAG T-3'

Ntrk1-R (16835): 5'-TCT ATG TGT GAG GTA TGT GCA TC-3'

**Reaction**

	Volume (μl)
H <sub>2</sub> O	2.5
Econo 2x mater mix	5
Ntrk1- F (10 μM)	1
Ntrk1- R(10 μM)	1
Template (Genomic DNA)	0.5
Final	10


**Use filter pipette tips when aliquoting reagents for PCR as well as when preparing master mix!**

**PCR program:**      95°C, 5 min  
                          95°C, 30 sec / 60°C, 30 sec / 72°C, 35 sec (30 cycles)  
                          72°C, 5 min

**Gel:** 2 % agarose gel, 135V, 30 min

**Results:** WT allele 472 bp; mutant allele ~600 bp

### **Genotyping of FGFR1 conditional KO mice**



Original paper: FGFR 1 is independently required in both developing mid-and hindbrain for sustained response to isthmic signals. The EMBO Journal, 2003, Vol. 22, 1811-1823.

#### **Primer:**

FGFR1-F: 5'- agggtcatcgaaatggacaag -3'

FGFR1-R: 5'- atcctcctgcttccttcaga -3'

#### **Reaction**

	Volume ( $\mu$ l)
H <sub>2</sub> O	2.5
Econo 2x mater mix	5
FGFR1- F (10 $\mu$ M)	1
FGFR1- R(10 $\mu$ M)	1
Template (Genomic DNA)	0.5
Final	10

**Use filter pipette tips when aliquoting reagents for PCR as well as when preparing master mix!**


**PCR program:**      95°C, 5 min  
                        95°C, 30 sec / 60°C, 30 sec / 72°C, 30 sec (34 cycles)  
                        72°C, 5 min  
                        \* the same program as FGF9

**Gel:** 2 % agarose gel, 135V, 30 min

**Results:** WT allel 200 bp; mutant allel 400 bp

### Genotyping of FGFR3 conditional KO mice

Original paper: Generation of Fgfr3 Conditional Knockout Mice. International Journal of Biological Sciences, 2010, 6(4):327-332.



#### **Primer:**

FGFR3-F: 5'- tgtaaagggtgggtggtag -3'

FGFR3-R: 5'- gctccctgtcctgcctcg -3'

#### **Reaction**

	Volume ( $\mu$ l)
H <sub>2</sub> O	2.5
Econo 2x mater mix	5
FGFR3- F (10 $\mu$ M)	1
FGFR3- R(10 $\mu$ M)	1
Template (Genomic DNA)	0.5
Final	10

**Use filter pipette tips when aliquoting reagents for PCR as well as when preparing master mix!**

**PCR program:** 95°C, 5 min

95°C, 30 sec / 60°C, 30 sec / 72°C, 30 sec (34 cycles)

72°C, 5 min

\* the same program as FGF9

**Gel:** 2 % agarose gel, 135V, 30 min

**Results:** WT allel 260 bp; mutant allel 320 bp

**Protocol-Genotyping**  
**Jessica Felker/Dr. Hui-Chen Lu's Lab**

**Genotyping of FGFR3-TDII (FGFR3 conditional) overexpression mice**

Fgfr3<sup>0</sup> with the K644E mutation, reflects the embryonic onset of neonatal lethal dwarfism, thanatophoric dysplasia type II (TDII)

Original paper: Human Molecular Genetics, 2000, Vol. 9, No. 11, 1603-1613

A neonatal lethal mutation in FGFR3 uncouples proliferation and differentiation of growth plate chondrocytes in embryos

**Primer:**

FGFR3-TDII-F: 5'- cgtggagttccactgcaag -3'

FGFR3-TDII -R: 5'- caccagccacgcagagtatg -3'

**Reaction**

	Volume (μl)
H <sub>2</sub> O	2.5
Econo 2x mater mix	5
FGFR3- F (10 μM)	1
FGFR3- R(10 μM)	1
Template (Genomic DNA)	0.5
Final	10

**Use filter pipette tips when aliquoting reagents for PCR as well as when preparing master mix!**

**PCR program:**      95°C, 5 min  
                          95°C, 30 sec / 58°C, 30 sec / 72°C, 25 sec (34 cycles)  
                          72°C, 5 min

**Gel:** 2 % agarose gel, 135V, 30 min


**Results:** mutant allel 250 bp

**Protocol-Genotyping**  
Jessica Felker/Dr. Hui-Chen Lu's Lab

**Genotyping of GAD67-GFP Knock-In mice**

Original paper: The Journal of Comparative Neurology, 2003, 467, 60-79.

Green Fluorescent Protein Expression and colocalization with calretinin, parvalbumin, and somatostatin in the GAD67-GFP Knock-in mouse



**Primer:**

GAD67-common (809-8): 5'-ggcacacagctccctctgtttgc -3'

GAD67-WT (809-9): 5'-gcttccttgcgttccgacag -3'

GAD67-KI (809-10): 5'-ctgtttcgccatgatatacg-3'

**Reaction**

	Volume (μl)
H <sub>2</sub> O	1.5
Econo 2x mater mix	5
GAD67-common (10 μM)	1
GAD67-WT (10 μM)	1
GAD67-KI (10 μM)	1
Template (Genomic DNA)	0.5
Final	10

**Use filter pipette tips when aliquoting reagents for PCR as well as when preparing master mix!**

**PCR program:** 95°C, 5 min


95°C, 30 sec / 60°C, 30 sec / 72°C, 25 sec (35 cycles)

72°C, 5 min

**Gel:** 2 % agarose gel, 135V, 30 min

**Results:** WT 265 bp; KI 564 bp

## **Genotyping of GPR55**



Wager-Miller, J., Wang, S., O'Leary, D., Monory, K., ... Lu, H.-C. (2010). Requirement of cannabinoid CB<sub>1</sub> receptor in cortical pyramidal neurons for appropriate development of corticothalamic and thalamocortical projections. *The European Journal of Neuroscience*, 32(5), 693–706.

<http://doi.org/10.1111/j.1460-9568.2010.07337.x>

**Primer:** GPR55-WT-F: 5'-GCCATCCAGTACCCGATCC-3'

GPR55-WT-R: 5'-GTCCAAGATAAAGCGGTTCC-3'

GPR55-Neo-F: 5'-GCAGCGCATCGCCTTCTATC-3'

GPR55-Neo-R: 5'-TCAAGCTACGTTTGGGTT-3'

Reactions:

GPR55-WT primer mix: GPR55-WT-F and GPR55-WT-R

GPR55-neo primer mix: GPR55-Neo-F and GPR55-Neo-R

**Reaction:**

	Volume ( $\mu$ l)
H <sub>2</sub> O	2
2X Econo Master Mix	5
Primer-F (10 $\mu$ M)	1
Primer-R (10 $\mu$ M)	1
Template (Genomic DNA)	1
Final	10

**Use filter pipette tips when aliquoting reagents for PCR as well as when preparing master mix!**

**PCR program:** 95°C, 5 min

95°C, 30 sec / 58°C, 30 sec / 72°C, 25 sec (35 cycles)

72°C, 5 min

**Gel:** 2% agarose gel, 135V, 30 min

**Results:** WT- 441 bp and Neo- 301 bp

**Protocol-Genotyping**  
Jessica Felker/Dr. Hui-Chen Lu's Lab

**Genotyping of Cx3cr1**

**Primer:** 12266 common: 5'-AAGACTCACGTGGACCTGCT-3'

14314 Mutant Reverse: 5'-CGGTTATTCAACTTGCACCA-3'

16221 Wild Type Reverse: 5'-AGGATGTTGACTTCCGAGTTG-3'

**Reaction:**

	Volume (μl)
H <sub>2</sub> O	1
2X Econo Master Mix	5
12266 common (10 μM)	1
14314 mutant reverse (10 μM)	1
16221 wild type reverse (10 μM)	1
Template (Genomic DNA)	1
Final	10

**Use filter pipette tips when aliquoting reagents for PCR as well as when preparing master mix!**

PCR program:	Cycling			
	Step #	Temp °C	Time	Note
1	94	2 min	-	
2	94	20sec	-	
3	65	15sec	-0.5 C per cycle	
4	68	10sec	-	
5	-	-	repeat steps 2-4 for 10 cycles	
6	94	15sec	-	
7	60	15sec	-	
8	72	10sec	-	
9	-	-	repeat steps 6-8 for 28 cycles	
10	72	2 min	-	
11	10	-	hold	

**Gel:** 2% agarose gel, 135V, 30 min

**Results:**

Mutant: ~300 bp

Heterozygous: ~300 bp and 695 bp

Wild type: 65 bp